


DISINFECTANT
PT2 BIOCIDES




UMONIUM³⁸

MASTER TISSUES

✓ FORMULA FREE FROM

- CMR* COMPONENT
- ENDOCRINE DISRUPTOR

SPORICIDE



“ Ultra-durable wet wipe with a broad-spectrum of activity ”

UMONIUM³⁸ **MASTER TISSUES** works within 5 seconds. Guaranteed performance for the daily treatment of rooms, shared rooms, equipment, as well as for *Clostridium difficile* situations (curative procedure).

UMONIUM³⁸ **MASTER TISSUES** is a biocidal product. Use biocides safely. Read the label and product information before use.

* Carcinogenic, Mutagenic or Reprotoxic

laboratoire
HUCKERT'S
INTERNATIONAL
mastered disinfection





UMONIUM³⁸ MASTER TISSUES

One 20 x 20cm wipe treats
an area of 3 m².



✓ **FORMULA FREE FROM**
→ CMR* COMPONENT
→ ENDOCRINE DISRUPTOR

FIELDS OF APPLICATION

- Furniture, high and low surfaces
- Health equipment
- Shoes, boots and goggles that have been in contact with contaminants

COMPOSITION

White non-woven wipes (20cm x 20cm) made of 50% wood pulp and 50% polyester and soaked in a solution of **UMONIUM³⁸ MASTER**: Alkyl (C12-16) dimethylbenzyl ammonium chloride (ADBAC/BKC (C12-16)) 11.9g/L.

REGULATORY STANDARDS

- The bactericidal, fungicidal and sporicidal activity of the soaked wipe has been validated according to standards ASTM 2967-15 and EN 16615
- The bactericidal, fungicidal, levoricidal, mycobactericidal and virucidal activity of the impregnation solution has been validated according to EN 13727, EN 14561, EN 13624, EN 14562, EN 14348, EN 14563 and EN 14476
- Quality and traceability guaranteed by ISO 9001:2015
- Complies with biocide regulation 528/2012

INSTRUCTIONS FOR USE

Cleaning/Pre-disinfection: using an evenly dampened wipe, rub the surface to remove any visual traces of organic material.


Disinfection: repeat the operation on a visually clean surface ensuring a minimum contact time of 5 seconds, only rinse with clean drinking water if necessary, leave it to dry.

PROPERTIES

- PT2 biocide
- Detergent- disinfectant wipe with a broad spectrum of activity
- Non-toxic: LD₅₀ (rat 5840 mg/kg), no dangerous vapours
- Non-corrosive: pH-neutral
- Ready-to-use product, non-irritating for the skin (no transdermal toxicity)
- Free from aldehyde, peroxide and biguanide, and other carcinogenic, mutagenic, or reprotoxic components
- Non-coagulant
- Non-flammable
- Bactericide, fungicide and sporicide activities
- 100% compatible on fibre optic, rubber, polycarbonate, acrylic compound, glass, Pyrex, PU, PVC, HDPE, PET, neoprene, latex, silicone, paints and stainless steel 410
- Effective, fast-acting and efficient (dissolves blood)
- Effective in the presence of organic materials
- Active within 5 seconds
- Residual on non-rinsed surfaces (bacteriostatic film)
- No residue after rinsing
- Stability: 36 months date of manufacture



PROCEDURE FOR USE



STEP 1: CLEANING/PRE-DISINFECTION



Take a wipe and proceed with wet wiping

If surfaces are visually clean, go to step 2

STEP 2: DISINFECTION



Take a new wipe and proceed with wet wiping



Contact time
min. 5 seconds



Rinse with clean
drinking water, if
necessary*

* For objects (such as toys) that may be put in the mouth by children

ITEM No.	NAME	BOX
PF 10750	Box of 100 wipes (20 x 20cm)	6 x 100

* Carcinogenic, Mutagenic or Reprotoxic

FP00412 EN - Ed. 3.2 - 03/2022 - Product photos are non-binding.



UMONIUM³⁸ is a registered trademark of Laboratoire Huckert's International